

Naif1 Rabbit Pab - Q69Yi7

Item Number: CM0018421

Introduction

Polyclonal antibody targeting the 35kDa nuclear apoptosis-inducing factor 1 (NAIF1) protein for WB and ELISA applications in human, mouse, and rat samples. Ideal for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	NAIF1 Rabbit pAb
Remarks/Alias	C9orf90; bA379C10.2; NAIF1
Species	Human
GeneID (Human)	203245
GeneID	203245
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human NAIF1 (NP_931045.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q69YI7
Protein Weight	35kDa
Shipping	Ice bag