


Nanog Rabbit Mab - Q80Z64

Item Number: CM0012146



Introduction

Rabbit monoclonal antibody against mouse Nanog (Homeobox protein NANOG). Validated for ELISA, WB, ChIP, and intracellular flow cytometry. Cross-reacts with human Nanog. Ideal for stem cell and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Nanog Rabbit mAb
Remarks/Alias	ENK; Stm1; ecat4; 2410002E02Rik
Species	Mouse
GeneID (Human)	79923
GeneID	71950
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 160-305 of mouse Nanog(NP_082292.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	ELISA, WB, ChIP, FC (intra)
Cross Reactivity	Human, Mouse
SWISS	Q80Z64
Protein Weight	34kDa
Shipping	Ice bag