

Napsa Rabbit Polyclonal Antibody For Wb, If-P, Elisa - O96009

Item Number: CM0026045

Introduction

High-quality NAPSA rabbit polyclonal antibody validated for WB, IF-P, and ELISA. Reacts with human, mouse, rat. Ideal for studying napsin-A function in pneumocyte surfactant processing.

[Learn More](#)

Parameter	Value
Product Name	NAPSA Rabbit pAb
Remarks/Alias	KAP; Kdap; NAP1; NAPA; SNAPA; NAPSA
Species	Human
GeneID	9476
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 321-420 of human NAPSA (NP_004842.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O96009
Protein Weight	45kDa
Shipping	Ice bag