


Nat10 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H0A0

Item Number: CM0026060



Introduction

High-quality NAT10 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Targets human and rat NAT10, a 116 kDa RNA cytidine acetyltransferase involved in mRNA and rRNA ac4C modification and ribosome biogenesis.

[Learn More](#)

Parameter	Value
Product Name	NAT10 Rabbit pAb
Remarks/Alias	ALP; Kre33; NET43; NAT10
Species	Human
Gene ID (Human)	55226
Gene ID	55226
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human NAT10 (NP_078938.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9H0A0
Protein Weight	116kDa
Shipping	Ice bag