

Nat2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P11245

Item Number: CM0013388



Introduction

NAT2 rabbit polyclonal antibody for WB, IHC-P, and ELISA applications. Reacts with human, mouse, rat NAT2 (SWISS: P11245). Suitable for drug metabolism, carcinogen bioactivation, and pharmacogenomic studies.

[Learn More](#)

Parameter	Value
Product Name	NAT2 Rabbit pAb
Remarks/Alias	AAC2; PNAT; NAT-2; NAT2
Species	Human
GeneID (human)	10
GeneID	10
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 171-290 of human NAT2 (NP_000006.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11245
Protein Weight	34kDa
Shipping	Ice bag