


Nav2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8IvL1

Item Number: CM0018806



Introduction

CamelBio offers NAV2 rabbit polyclonal antibody (SWISS Q8IVL1) specific for human NAV2.

Validated for WB, IHC-P, ELISA; detects 268kDa protein involved in neuronal development and cytoskeleton organization.

[Learn More](#)

Parameter	Value
Product Name	NAV2 Rabbit pAb
Remarks/Alias	HELAD1; RAINB1; POMFIL2; UNC53H2; STEERIN2; NAV2
Species	Human
Gene ID (Human)	89797
Gene ID	89797
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1800-1980 of human NAV2 (NP_660093.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q8IVL1
Protein Weight	268kDa
Shipping	Ice bag