


Ncoa4 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q13772

Item Number: CM0004523



Introduction

NCOA4 Rabbit Monoclonal Antibody (Q13772) validated for WB, IP, ELISA with cross-reactivity in Human, Mouse, Rat. Ideal for ferritinophagy, iron homeostasis, and nuclear receptor coactivator studies.

[Learn More](#)

Parameter	Value
Product Name	NCOA4 Rabbit mAb
Remarks/Alias	RFG; ELE1; PTC3; ARA70; NCOA4
Species	Human
Gene ID (Human)	8031
Gene ID	8031
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13772
Protein Weight	70kDa
Shipping	Ice bag