


Ncr3 Rabbit Pab For Wb, Elisa - O14931

Item Number: CM0029910



Introduction

High-quality NCR3 Rabbit Polyclonal Antibody (pAb) for WB, ELISA. Reactive with human NCR3 (NKp30, 22kDa). Ideal for NK cell cytotoxicity and cancer immunology research.

[Learn More](#)

Parameter	Value
Product Name	NCR3 Rabbit pAb
Remarks/Alias	1C7; MALS; CD337; LY117; NKp30; NCR3
Species	Human
GeneID (human)	259197
GeneID	259197
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-135 of human NCR3 (NP_667341.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14931
Protein Weight	22kDa
Shipping	Ice bag