


Ncs1 Rabbit Monoclonal Antibody For Wb And Elisa - P62166

Item Number: CM0002129



Introduction

Recombinant NCS1 rabbit monoclonal antibody validated for Western blot and ELISA. Detects human, mouse, and rat neuronal calcium sensor 1 with a predicted 22 kDa band. Ideal for neuroscience and calcium signaling research.

[Learn More](#)

Parameter	Value
Product Name	NCS1 Rabbit mAb
Remarks/Alias	FLUP; FREQ; NCS1
Species	Human
GeneID (Human)	23413
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 91-190 of human NCS1 (P62166).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62166
Protein Weight	22kDa
Shipping	Ice bag