

Ndp Rabbit Pab - Q00604

Item Number: CM0020859



Introduction

Rabbit polyclonal antibody against Norrin (NDP) for Western blot and ELISA applications. Reacts with human NDP. Ideal for retinal vascularization and Wnt signaling research.

[Learn More](#)

Parameter	Value
Product Name	NDP Rabbit pAb
Remarks/Alias	ND; EVR2; FEVR; NDP
Species	Human
GeneID (Human)	4693
GeneID	4693
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-133 of human NDP (NP_000257.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q00604
Protein Weight	15kDa
Shipping	Ice bag