

Ndufa4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00483

Item Number: CM0020849



Introduction

High-quality NDUFA4 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat samples. Ideal for mitochondrial respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	NDUFA4 Rabbit pAb
Remarks/Alias	MLRQ; CI-9k; COXFA4; MISTR1; MRCAF1; CI-MLRQ; MC4DN21; NDUFA4
Species	Human
GeneID (Human)	4697
GeneID	4697
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-81 of human NDUFA4 (NP_002480.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00483
Protein Weight	9kDa
Shipping	Ice bag