


Ndufb2 Rabbit Pab For Wb, Elisa - O95178

Item Number: CM0028220



Introduction

Rabbit polyclonal antibody targeting human NDUFB2, a Complex I accessory subunit. Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for mitochondrial respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	NDUFB2 Rabbit pAb
Remarks/Alias	AGGG; CI-AGGG; NDUFB2
Species	Human
GeneID (human)	4708
GeneID	4708
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-105 of human NDUFB2 (NP_004537.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95178
Protein Weight	12kDa
Shipping	Ice bag