


Ndufb5 Rabbit Polyclonal Antibody For Wb, Elisa - O43674

Item Number: CM0029429



Introduction

NDUFB5 Rabbit pAb detects human NDUFB5, an accessory subunit of mitochondrial Complex I, with cross-reactivity to mouse and rat. Validated for WB and ELISA. Ideal for respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	NDUFB5 Rabbit pAb
Remarks/Alias	SGDH; CISGDH; NDUFB5
Species	Human
Gene ID (Human)	4711
Gene ID	4711
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-189 of human NDUFB5 (NP_002483.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43674
Protein Weight	22kDa
Shipping	Ice bag