


Ndufs4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43181

Item Number: CM0025220



Introduction

NDUF54 Rabbit Polyclonal Antibody validated for WB, IHC-P, and ELISA. Targets the 20 kDa mitochondrial Complex I subunit. Cross-reacts with human, mouse, and rat. Ideal for respiratory chain and mitochondrial disease studies.

[Learn More](#)

Parameter	Specification
Product Name	NDUF54 Rabbit pAb
Remarks/Alias	AQDQ; CI-18; MC1DN1; CI-AQDQ; CI-18 kDa; NDUF54
Species	Human
GeneID (Human)	4724
GeneID	4724
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human NDUF54 (NP_002486.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43181
Protein Weight	20kDa
Shipping	Ice bag