

Neurl1B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - A8Mq27

Item Number: CM0019866

Introduction

NEURL1B rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets E3 ubiquitin-protein ligase NEURL1B, a Notch pathway regulator with cytoplasmic localization.

[Learn More](#)

Parameter	Value
Product Name	NEURL1B Rabbit pAb
Remarks/Alias	neur2; NEURL3; RNF67B; hNeur2; NEURL1B
Species	Human
GeneID (Human)	54492
GeneID	54492
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-440 of human NEURL1B (NP_001136123.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	A8MQ27
Protein Weight	59kDa
Shipping	Ice bag