

Neurod1 Rabbit Monoclonal Antibody For Wb And Elisa - Q13562

Item Number: CM0003005

Introduction

High-quality rabbit monoclonal antibody against human NeuroD1 (bHLHa3, BETA2). Validated for Western blot and ELISA, with cross-reactivity to mouse. Suitable for neurogenesis and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	NeuroD1 Rabbit mAb
Remarks/Alias	T2D; BETA2; BHF-1; MODY6; NEUROD; bHLHa3; NeuroD1
Species	Human
GeneID (Human)	4760
GeneID	4760
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NeuroD1 (NP_002491.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q13562
Protein Weight	40kDa
Shipping	Ice bag