

Neurofilament H Rabbit Monoclonal Antibody For Wb - P12036

Item Number: CM0006828

Introduction

Rabbit monoclonal antibody against human Neurofilament H (NFH), validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for neuronal research and neurodegenerative disease studies.

[Learn More](#)

Parameter	Value
Product Name	Neurofilament H Rabbit mAb
Remarks/Alias	NFH; CMT2CC; Neurofilament H
Species	Human
Gene ID	4744
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 873-1020 of human Neurofilament H. (NP_066554.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12036
Protein Weight	112kDa
Shipping	Ice bag