

Neurofilament M Rabbit Mab For Wb, If/Icc, Elisa - P07197

Item Number: CM0002712

Introduction

Rabbit monoclonal antibody against human Neurofilament M (NFM, NEF3, Neurofilament 3). Validated for Western blot, immunofluorescence/ICC, and ELISA. Ideal for studying neuronal cytoskeleton and neurofilament organization.

[Learn More](#)

Parameter	Value
Product Name	Neurofilament M Rabbit mAb
Remarks/Alias	NFM; NEF3; NF-M; Neurofilament M
Species	Human
GeneID (Human)	4741
GeneID	4741
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 600-700 of human Neurofilament M (P07197).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P07197
Protein Weight	102 kDa
Shipping	Ice bag