

Neurog1 Rabbit Pab - Q92886

Item Number: CM0025902

Introduction

NEUROG1 rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human and mouse. Targets a 26kDa transcriptional regulator involved in neuronal differentiation. Research use only.

[Learn More](#)

Parameter	Value
Product Name	NEUROG1 Rabbit pAb
Remarks/Alias	AKA; ngn1; Math4C; bHLHa6; NEUROD3; NEUROG1
Species	Human
GeneID (Human)	4762
GeneID	4762
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-237 of human NEUROG1 (Q92886)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92886
Protein Weight	26kDa
Shipping	Ice bag