


Nf-Kb P65/Rela Mouse Monoclonal Antibody For Wb, Elisa - Q04206

Item Number: CM0003843



Introduction

High-quality mouse monoclonal antibody against human NF-kB p65/RelA (SWISS Q04206) validated for Western blot and ELISA. Suitable for human, mouse, and rat samples. Ideal for inflammation and immunology research.

[Learn More](#)

Parameter	Value
Product Name	NF-kB p65/RelA Mouse mAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; NF-kB p65/RelA
Species	Human
GeneID (Human)	5970
GeneID	5970
Immunogen	Recombinant protein of human NF-kB p65
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04206
Protein Weight	60 kDa
Shipping	Ice bag