


Nf-Kb P65/Rela Rabbit Polyclonal Antibody - Q04206

Item Number: CM0020593



Introduction

Rabbit polyclonal antibody targeting the p65/RelA subunit of NF-kB, validated for WB, IF/ICC, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for inflammation, immunity, and transcriptional regulation studies. - Q04206

[Learn More](#)

Parameter	Value
Product Name	NF-kB p65/RelA Rabbit pAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; NF-kB p65/RelA
Species	Human
GeneID (Human)	5970
GeneID	5970
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-550 of human NF-kB p65/RelA (NP_068810.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04206
Protein Weight	60 kDa
Shipping	Ice bag