


Nf-Kb2 Rabbit Pab - Polyclonal Antibody For Wb, If/Icc, Elisa - Q00653

Item Number: CM0015368



Introduction

NF-κB2 Rabbit Polyclonal Antibody (pAb) detects p100/p52 subunits in WB, IF/ICC, ELISA. Reacts with human and mouse. Suitable for NF-κB pathway research, immune response, and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	NF-κB2 Rabbit pAb
Remarks/Alias	p52; p100; H2TF1; LYT10; CVID10; LYT-10; NF-κB2; p49/p100; NF-κB2
Species	Human
GeneID (Human)	4791
GeneID	4791
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 690-899 of human NF-κB2 (NP_002493.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q00653
Protein Weight	97kDa
Shipping	Ice bag