


Nf-Kb2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q00653

Item Number: CM0026282



Introduction

Rabbit polyclonal antibody targeting human NF-κB2/p100/p52. Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying NF-kappa-B signaling, inflammation, and immunity.

[Learn More](#)

Parameter	Specification
Product Name	NF-κB2 Rabbit pAb
Remarks/Alias	p52; p100; H2TF1; LYT10; CVID10; LYT-10; NF-κB2; p49/p100; NF-κB2
Species	Human
Gene ID (Human)	4791
Gene ID	4791
Immunogen	Recombinant protein corresponding to amino acids 690-899 of human NF-κB2 (NP_002493.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00653
Protein Weight	97kDa
Shipping	Ice bag