


Nfasc Rabbit Pab - O94856

Item Number: CM0030834



Introduction

Rabbit polyclonal antibody targeting human Neurofascin (NFASC), a 150kDa cell adhesion protein. Validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for neuroscience research on neurite extension, axonal guidance, and myelination.

[Learn More](#)

Parameter	Value
Product Name	NFASC Rabbit pAb
Remarks/Alias	NF; NRCAML; NEDCPMD; NFASC
Species	Human
GeneID (human)	23114
GeneID	23114
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94856
Protein Weight	150kDa
Shipping	Ice bag