


Nfatc2 Rabbit Pab - Q13469

Item Number: CM0030417



Introduction

NFATC2 Rabbit pAb is a polyclonal antibody raised against recombinant human NFATC2 protein, suitable for WB, IF/ICC, IF-P, IHC-P, ELISA. It detects NFATC2 (100 kDa) in human, mouse, rat samples and is ideal for studying NFAT signaling in T-cell activation.

[Learn More](#)

Parameter	Value
Product Name	NFATC2 Rabbit pAb
Remarks/Alias	JCOSL; NFAT1; NFATP; NFATC2
Species	Human
GeneID (Human)	4773
GeneID	4773
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human NFATC2 (NP_775114.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13469
Protein Weight	100kDa
Shipping	Ice bag