


Nfatc4 Rabbit Pab For Wb And Elisa - Q14934

Item Number: CM0024542



Introduction

Rabbit polyclonal antibody targeting NFATC4 (NFAT3), validated for Western blot and ELISA. Detects human and mouse NFATC4, a transcription factor involved in immune, cardiac, and neural development.

[Learn More](#)

Parameter	Value
Product Name	NFATC4 Rabbit pAb
Remarks/Alias	NFAT3; NF-AT3; NF-ATC4; NFATC4
Species	Human
GeneID (Human)	4776
GeneID	4776
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 560-710 of human NFATC4 (NP_004545.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14934
Protein Weight	95kDa
Shipping	Ice bag