


Nfkb1 Rabbit Pab For Wb, If/Icc, Elisa - P19838

Item Number: CM0015082



Introduction

NFKB1 Rabbit pAb (p105/p50) for WB, IF/ICC, ELISA. Detects human, mouse, rat NFKB1, a key transcription factor in inflammation and immunity. Ideal for signaling and disease research.

[Learn More](#)

Parameter	Value
Product Name	NFKB1 Rabbit pAb
Remarks/Alias	KBF1; EBP-1; NF-kB; CVID12; NF-kB1; NFKB-p50; NFkappaB; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta; NFKB1
Species	Human
GeneID (human)	4790
GeneID	4790
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 41-365 of human NFKB1 (NP_001158884.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19838
Protein Weight	105kDa
Shipping	Ice bag