

Nfkb1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P19838

Item Number: CM0026003



Introduction

NFKB1 Rabbit Polyclonal Antibody (P19838) validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Targets NF-kB p105 precursor and p50 subunit, key regulators of immune and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	NFKB1 Rabbit pAb
Remarks/Alias	KBF1; EBP-1; NF-kB; CVID12; NF-kB1; NFKB-p50; NFkappaB; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta; NFKB1
Species	Human
GeneID (Human)	4790
GeneID	4790
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 41-365 of human NFKB1 (NP_001158884.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19838
Protein Weight	105kDa
Shipping	Ice bag