

Nfya Rabbit Mab For Wb, Ip, Ihc-P, Chip, Elisa - P23511

Item Number: CM0011041



Introduction

Recombinant rabbit monoclonal antibody against NFYA (NF-YA), validated for WB, IP, IHC-P, ChIP, and ELISA. Reacts with human, mouse, rat. Ideal for transcription and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	NFYA Rabbit mAb
Remarks/Alias	HAP2; CBF-A; CBF-B; NF-YA
Species	Human
Gene ID (Human)	4800
Gene ID	4800
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 262-332 of human NFYA (NP_002496.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23511
Protein Weight	34kDa/37kDa
Shipping	Ice bag