

Ngef Rabbit Polyclonal Antibody For Wb, If/Icc - Q8N5V2

Item Number: CM0021217

Introduction

NGEF Rabbit Polyclonal Antibody validated for Western blot, immunofluorescence/ICC, and ELISA. Detects endogenous levels of Ephexin-1 (NGEF) in human, mouse, and rat samples. Ideal for axon guidance and neurogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	NGEF Rabbit pAb
Remarks/Alias	EPHEXIN; ARHGEF27; NGEF
Species	Human
Gene ID (Human)	25791
Gene ID	25791
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 429-710 of human NGEF (NP_062824.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N5V2
Protein Weight	82kDa
Shipping	Ice bag