

Nlrp1 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - Q9C000

Item Number: CM0019951



Introduction

Rabbit polyclonal antibody targeting human NLRP1 (SWISS: Q9C000), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for innate immunity, inflammasome, and pyroptosis studies.

[Learn More](#)

Parameter	Value
Product Name	NLRP1 Rabbit pAb
Remarks/Alias	NAC; JRRP; MSPC; AIADK; CARD7; CIDED; NALP1; SLEV1; DEFCAP; PP1044; VAMAS1; CLR17.1; DEFCAP-L/S; NLRP1
Species	Human
GeneID (Human)	22861
GeneID	22861
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1370-1473 of human NLRP1 (NP_127497.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C000
Protein Weight	166kDa
Shipping	Ice bag