


Nlrp3 Rabbit Mab - Q8R4B8

Item Number: CM0007662



Introduction

Anti-NLRP3 Rabbit monoclonal antibody for mouse samples, suitable for WB, IF/ICC, IP, and ELISA. NLRP3 is a key inflammasome sensor mediating IL-1 β secretion and pyroptosis. Ideal for innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	NLRP3 Rabbit mAb
Remarks/Alias	FCU; MWS; FCAS; Cias1; Mmig1; NALP3; Pypaf1; Ail/AVP; AGTAVPRL; NLRP3
Species	Mouse
GeneID (Human)	114548
GeneID	216799
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of mouse NLRP3(NP_665826.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Mouse
SWISS	Q8R4B8
Protein Weight	118kDa
Shipping	Ice bag