


Nlrp3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8R4B8

Item Number: CM0022503



Introduction

NLRP3 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Raised against mouse Nlrp3 aa 1-93. Explore NLRP3 inflammasome pathways with this 118kDa target.

[Learn More](#)

Parameter	Value
Product Name	NLRP3 Rabbit pAb
Remarks/Alias	FCU; MWS; FCAS; Cias1; Mmig1; NALP3; Pypaf1; Ail/AVP; AGTAVPRL; NLRP3
Species	Mouse
GeneID (human)	216799
GeneID	216799
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of mouse NLRP3 (NP_665826.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8R4B8
Protein Weight	118kDa
Shipping	Ice bag