


Nlrp3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96P20

Item Number: CM0026163



Introduction

High-quality rabbit polyclonal NLRP3 antibody for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat NLRP3 (~118 kDa). Ideal for studying inflammasome activation, IL-1 β secretion, and autoinflammatory disorders.

[Learn More](#)

Parameter	Value
Product Name	NLRP3 Rabbit pAb
Remarks/Alias	AIL; AVP; FCU; MWS; FCAS; KEFH; CIAS1; FCAS1; NALP3; C1orf7; CLR1.1; DFNA34; PYPAF1; AGTAVPRL; NLRP3
Species	Human
Gene ID	114548
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-210 of human NLRP3 (NP_004886.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96P20
Protein Weight	118 kDa
Shipping	Ice bag