


NlrX1 Rabbit Pab - Q86Ut6

Item Number: CM0026586



Introduction

NLRX1 Rabbit pAb is a polyclonal antibody validated for WB, IP, ELISA. It detects a 108kDa mitochondrial protein in human, mouse, rat, and serves as a negative regulator of antiviral innate immunity.

[Learn More](#)

Parameter	Value
Product Name	NLRX1 Rabbit pAb
Remarks/Alias	NOD5; NOD9; NOD26; DLNB26; CLR11.3; NLRX1
Species	Human
Gene ID (Human)	79671
Gene ID	79671
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 87-360 of human NLRX1 (NP_078894.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86UT6
Protein Weight	108kDa
Shipping	Ice bag