


Nm23 Rabbit Polyclonal Antibody For IHC-P And Elisa - P15531 / P22392

Item Number: CM0023973



Introduction

Rabbit polyclonal antibody targeting human nm23 (NME1/NME2) for IHC-P and ELISA applications. Detects human and rat nm23. Suitable for nucleoside diphosphate kinase, metastasis, and cell motility research.

[Learn More](#)

Parameter	Value
Product Name	nm23 Rabbit pAb
Remarks/Alias	NB; AWD; NBS; GAAD; NDKA; NM23; NDPKA; NDPK-A; NM23-H1; PUF; NDKB; NDPKB; NM23B; NDPK-B; NM23-H2; nm23
Species	Human
Gene ID	4830/4831
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-152 of human nm23 (NP_000260.1/NP_002503.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P15531/P22392
Protein Weight	17kDa
Shipping	Ice bag