


Nme4 Rabbit Pab - 000746

Item Number: CM0028121



Introduction

NME4 Rabbit pAb is a polyclonal antibody against human NME4, suitable for WB and ELISA. It cross-reacts with human, mouse, and rat, and is valuable for studying mitochondrial nucleoside diphosphate kinase and phospholipid transfer functions.

[Learn More](#)

Parameter	Value
Product Name	NME4 Rabbit pAb
Remarks/Alias	NDPK-D; NM23H4; nm23-H4; NME4
Species	Human
Gene ID (Human)	4833
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-187 of human NME4 (NP_005000.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00746
Protein Weight	21kDa
Shipping	Ice bag