


Non-Phospho-Yap1-S127 Rabbit Pab For Wb, Ihc-P, Elisa - P46937

Item Number: CM0020946



Introduction

Rabbit polyclonal antibody recognizing non-phosphorylated YAP1 at Ser127; validated for WB, IHC-P, and ELISA in human, mouse, and rat. Essential for Hippo signaling and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	Non-phospho-YAP1-S127 Rabbit pAb
Remarks/Alias	YAP; YKI; COB1; YAP2; YAP-1; YAP65; Non-phospho-YAP1-S127
Species	Human
GeneID (Human)	10413
GeneID	10413
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human YAP1 (NP_001123617.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46937
Protein Weight	54kDa
Shipping	Ice bag