


## Nop10 Rabbit Pab - Q9Npe3

Item Number: CM0021660



### Introduction

Rabbit polyclonal antibody against NOP10 (H/ACA ribonucleoprotein subunit 3). Validated for WB and ELISA, reacts with human and mouse. For ribosome biogenesis and telomere maintenance research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | NOP10 Rabbit pAb                                                                                                |
| Remarks/Alias    | DKCB1; NOLA3; NOP10P; NOP10                                                                                     |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 55505                                                                                                           |
| GeneID           | 55505                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-64 of human NOP10 (NP_061118.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, ELISA                                                                                                       |
| Cross Reactivity | Human, Mouse                                                                                                    |
| SWISS            | Q9NPE3                                                                                                          |
| Protein Weight   | 8kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                         |