


Nop58 Rabbit Pab For Wb, If/Icc, Elisa - Q9Y2X3

Item Number: CM0027576



Introduction

High specificity NOP58 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA across human, mouse, and rat samples; ideal for ribosome biogenesis and snoRNP research.

[Learn More](#)

Parameter	Value
Product Name	NOP58 Rabbit pAb
Remarks/Alias	NOP5; HSPC120; NOP5/NOP58; NOP58
Species	Human
Gene ID (Human)	51602
Gene ID	51602
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 360-529 of human NOP58 (NP_057018.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2X3
Protein Weight	60kDa
Shipping	Ice bag