

Notch3 Rabbit Pab - Q9Um47

Item Number: CM0026607



Introduction

Rabbit polyclonal antibody against NOTCH3 for WB, IF/ICC, and ELISA applications. Validated in human, mouse, rat. Ideal for investigating Notch signaling, cell-fate determination, and CADASIL disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	NOTCH3 Rabbit pAb
Remarks/Alias	IMF2; LMNS; CASIL; CADASIL; CADASIL1; NOTCH3
Species	Human
GeneID (human)	4854
GeneID	4854
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1520-1640 of human NOTCH3 (NP_000426.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UM47
Protein Weight	244kDa
Shipping	Ice bag