

Nox5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Ph1

Item Number: CM0026535



Introduction

Rabbit polyclonal antibody targeting human NOX5, a calcium-dependent NADPH oxidase. Validated for WB, IF/ICC, and ELISA. Ideal for studying superoxide generation, ROS signaling, and angiogenesis in endothelial and immune cells.

[Learn More](#)

Parameter	Value
Product Name	NOX5 Rabbit pAb
Remarks/Alias	NOX5
Species	Human
GeneID (Human)	79400
GeneID	79400
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 661-765 of human NOX5 (NP_078781.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q96PH1
Protein Weight	86kDa
Shipping	Ice bag