


Noxa1 Rabbit Pab For Wb And Elisa - Q86Ur1

Item Number: CM0014030



Introduction

NOXA1 Rabbit pAb is a polyclonal antibody validated for Western blotting and ELISA, cross-reacting with human, mouse, and rat. It targets NOXA1, a key NADPH oxidase activator in ROS production, supporting studies in oxidative stress, host defense, and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	NOXA1 Rabbit pAb
Remarks/Alias	p51NOX; NY-CO-31; SDCCAG31; NOXA1
Species	Human
Gene ID (Human)	10811
Gene ID	10811
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86UR1
Protein Weight	51kDa
Shipping	Ice bag