


Nphp3 Rabbit Polyclonal Antibody For Wb And Elisa - Q7Z494

Item Number: CM0025379



Introduction

Rabbit polyclonal antibody against human NPHP3 (Nephrocystin-3), validated for WB and ELISA. Recombinant immunogen (aa1-130). Cross-reacts with mouse. For ciliopathy and Wnt signaling research.

[Learn More](#)

Parameter	Value
Product Name	NPHP3 Rabbit pAb
Remarks/Alias	MKS7; NPH3; RHPD; RHPD1; SLSN3; CFAP31; NPHP3
Species	Human
GeneID (human)	27031
GeneID	27031
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human NPHP3 (NP_694972.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z494
Protein Weight	151kDa
Shipping	Ice bag