

Npm1 (B23, Npm) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P06748

Item Number: CM0020671



Introduction

Rabbit polyclonal antibody targeting human NPM1 (Nucleophosmin, B23). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. For research on ribosome biogenesis, centrosome duplication, and cancer.

[Learn More](#)

Parameter	Value
Product Name	NPM1 Rabbit pAb
Remarks/Alias	B23, NPM, NPM1
Species	Human
GeneID (human)	4869
GeneID	4869
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 9-122 of human NPM1 (NP_002511.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06748
Protein Weight	33 kDa
Shipping	Ice bag