


Nppa Rabbit Pab - P01160

Item Number: CM0023410



Introduction

NPPA Rabbit pAb is a polyclonal antibody validated for WB, IF-P, and ELISA. It recognizes human NPPA (ANP) and cross-reacts with mouse and rat, ideal for cardiovascular and renal research.

[Learn More](#)

Parameter	Value
Product Name	NPPA Rabbit pAb
Remarks/Alias	ANF; ANP; CDD; CDP; PND; ATFB6; ATRST2; CDD-ANF; NPPA
Species	Human
GeneID (Human)	4878
GeneID	4878
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-151 of human NPPA (NP_006154.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01160
Protein Weight	16 kDa
Shipping	Ice bag