


Npr3 Rabbit Pab For Wb, If/Icc, Elisa - P17342

Item Number: CM0024259



Introduction

Rabbit polyclonal antibody against human NPR3 (natriuretic peptide clearance receptor). Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Suitable for studying cardiovascular and bone biology signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	NPR3 Rabbit pAb
Remarks/Alias	NPRC; ANP-C; ANPRC; BOMOS; NPR-C; ANPR-C; GUCY2B; C5orf23; NPR3
Species	Human
Gene ID (Human)	4883
Gene ID	4883
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-480 of human NPR3 (NP_001191304.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17342
Protein Weight	60kDa
Shipping	Ice bag