

Npy Rabbit Mab For Western Blot, Immunofluorescence, Immunohistochemistry, Elisa - P01303

Item Number: CM0004475



Introduction

Anti-NPY rabbit monoclonal antibody for WB, IF, IHC, ELISA. Detects human, mouse, and rat Neuropeptide Y (PYY4). Perfect for neurobiology and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	NPY Rabbit mAb
Remarks/Alias	PYY4
Species	Human
GeneID (Human)	4852
GeneID	4852
Immunogen	A synthesized peptide derived from human NPY
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01303
Protein Weight	11 kDa
Shipping	Ice bag