

Nrn1L Rabbit Pab For Wb, Elisa - Q496H8

Item Number: CM0021877

Introduction

NRN1L Rabbit Polyclonal Antibody (pAb) for WB and ELISA, validated in human, mouse, rat. 18kDa neuritin-like protein antibody for cell membrane and GPI-anchor research.

[Learn More](#)

Parameter	Value
Product Name	NRN1L Rabbit pAb
Remarks/Alias	UNQ2446; cpg15-2; MRCC2446; NRN1L
Species	Human
Gene ID (Human)	123904
Gene ID	123904
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-139 of human NRN1L (Q496H8)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q496H8
Protein Weight	18kDa
Shipping	Ice bag