


Ntan1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Ab6

Item Number: CM0013366



Introduction

Rabbit polyclonal antibody targeting human NTAN1 (PNAA/PNAD), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying N-end rule protein degradation.

[Learn More](#)

Parameter	Value
Product Name	NTAN1 Rabbit pAb
Remarks/Alias	PNAA; PNAD; NTAN1
Species	Human
GeneID (Human)	123803
GeneID	123803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human NTAN1 (NP_775745.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AB6
Protein Weight	35kDa
Shipping	Ice bag